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The novel, peptide-based hemostat

cardiovascular case report | Vol. 2

Haemostasis using PuraStat® during Left Ventricular Assist Device Implantation (LVAD)

CASE PRESENTATION

DIAGNOSIS
Dilated Cardiomyopathy, requiring Left Ventricular Assist Device implantation as a bridge to cardiac transplant.

PATIENT DEMOGRAPHICS
Female, 42 years

PROCEDURE

- Implantation of Left Ventricular Assist Device via thoracotomy approach
- Sewing ring attached to apex of left ventricle (Fig. 1)
- Outflow graft attached to ascending aorta
- PuraStat applied to apical sewing ring and anastomosis of outflow graft to aorta (Fig. 2 and Fig. 3)

POSTOPERATIVE COURSE
Uneventful
No postoperative bleeding complications
Extubated within 72 hours

FEEDBACK ON PURASTAT USAGE

André R. Simon, M.D., Ph.D.
Director of the Department for Heart and Lung Transplantation and Mechanical Circulatory Support.
Harefield Hospital, United Kingdom

"PuraStat did what it was supposed to do."

"I've already used PuraStat for over 60 cardiothoracic cases and PuraStat is now a standard of care for most of my cases. PuraStat can be used safely in Ventricular Assist Device implantation. It is easily applied around the sewing ring and the aorta to quickly achieve haemostasis. Even after application of PuraStat, it doesn't become hard, which allows you to place a suture right through the material when you need additional stitches to the site."

"The efficacy of PuraStat compares favourably with other haemostatic agents (glues, flowables) and the benefits of PuraStat such as transparency, ease of use, safety and better control are indeed satisfactory."

CR CV 002 GB V1 2016 09 20 Page 2 of 2

PuraStat®

The novel, peptide-based hemostat

cardiovascular case report | Vol. 3

Haemostasis using PuraStat® during 'Dor procedure'

CASE PRESENTATION

DIAGNOSIS
Giant left ventricular aneurysm after myocardial infarction (antero-septo-apical region), left ventricular ejection fraction 37%

PATIENT DEMOGRAPHICS
Male, 55 years

PROCEDURE
'Dor Procedure':

- Cannulation: Ascending aorta and Right atrium / inferior vena cava, full heparinization (300 IU/kg)
- Pump time: 68 min
- Aortic cross-clamp time: 37 min
- Cold blood cardioplegia ('Buckberg cardioplegia')
- Endoventricular patch used to complete the repair
- PuraStat was applied to the reconstructed left ventricular wall and on the vulnerable surface of the left ventricle

POSTOPERATIVE COURSE
Uneventful

FEEDBACK ON PURASTAT USAGE

Igor Rudez, M.D., Ph.D.
Assistant Professor
Head, Division of Aortic and Coronary Surgery Department of Cardiac and Transplant Surgery
University Hospital Dubrava, Zagreb, Croatia

"PuraStat is an excellent haemostatic tool, particularly in management of large, vulnerable, oozing surfaces. Ready-to-use off-the-shelf, no need for "cooking", warming, or mixing components. Completely transparent it does not blur the bleeding site."

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The novel, peptide-based hemostat

cardiovascular case report | Vol. 5

Haemostasis using PuraStat® in total aortic arch replacement

CASE PRESENTATION

DIAGNOSIS

Coronary artery disease and acute aortic dissection (type B)

PATIENT DEMOGRAPHICS

Male, 64 years

PROCEDURE

- 'Frozen Elephant Trunk procedure' with separated re-implantation of epi-aortic branches
- Double coronary artery bypass grafting (CABG) was performed
- PuraStat has been applied in 3 different areas during the procedure:
 - Distal anastomosis between polyester graft and descending aorta (Fig. 1)
 - Left sub-clavian artery anastomosis
 - Anastomosis between polyester graft and ascending aorta (Fig. 2)

POSTOPERATIVE COURSE

Uneventful.

Patient was weaned from cardiopulmonary bypass. No bleeding occurred.

Fig. 1: PuraStat was applied on distal anastomosis between polyester graft and descending aorta

Fig. 2: PuraStat was applied on anastomosis between polyester graft and ascending aorta

Fig. 3: Haemostasis achieved

FEEDBACK ON PURASTAT USAGE

Roberto Di Bartolomeo, M.D., Professor
Professor of Cardiac Surgery at Alma Mater Studiorum Bologna University
Director of Cardiac Surgery Unit
President of the Italian Society of Cardiac Surgery (SICCH)
Policlinico di S. Orsola - Azienda Ospedaliero-Universitaria di Bologna, Italy

"PuraStat is a safe and very effective hemostat device. It is very easy to prepare and to manage during the surgical operation. It is also very fast and it is able to stop bleeding in a few seconds. Its transparency helps a lot in case of secondary bleedings. In this particular case, PuraStat was applied in aorta anastomosis areas stopping any bleeding between the natural vessel and the prosthesis. I have also appreciated the fact that it can be used in different times during the operation without the need of changing the applicator."

CR CV 005 GB V1 2016 09 20 Page 2 of 2

PuraStat®

The novel, peptide-based hemostat

cardiovascular case report | Vol. 4

Haemostasis using PuraStat® in repair of ascending aorta aneurism

CASE PRESENTATION

DIAGNOSIS

Ascending aorta aneurism without aortic insufficiency

PATIENT DEMOGRAPHICS

Male, 56 years

History of arterial hypertension and obesity.

PROCEDURE

Using cardiopulmonary bypass and cardiac arrest:

- Replacement of ascending aorta with aortic valve conservation and coronary arteries re-implantation ('David procedure')
- The ascending aorta was replaced using a 30 mm polyester graft (Fig. 2)
- Aortic valve was then re-implanted inside the graft with continuous sutures
- The next step was coronary artery re-implantation. Finally the distal anastomosis with distal aorta was performed
- One application of PuraStat was used on prosthesis needle holes (Fig. 3)
- Another application of PuraStat was performed on the sutures before aortic declamping (Fig. 4)
- Good haemostasis was observed and approximately 2 ml of PuraStat was used in total

POSTOPERATIVE COURSE

Uneventful

No bleeding was observed in the intensive care unit and the patient did not need blood transfusion

Fig. 1: Resection of aneurism, preparation of coronary ostia and placement of stitches under aortic annulus


Fig. 2: The ascending aorta was replaced using a 30 mm polyester graft

Fig. 3: PuraStat was used on prosthesis needle holes

Fig. 4: Another application of PuraStat was performed on the sutures before aortic declamping

FEEDBACK ON PURASTAT USAGE

Fabien Doguet, M.D., Ph.D., Professor
Department of Thoracic and Cardiovascular Surgery

CHU-Hôpitaux de Rouen, France

"PuraStat is really easy to use both for the surgeon and for the nursing staff, who prepared the syringe. In ascending aorta surgery, where we need to do lots of sutures, the risk of bleeding is high, and the use of PuraStat is helpful to reduce operative bleeding"

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3-D Matrix, Ltd